


AirClean® Systems

"THE FUME CONTROL EXPERTS"

Performance & Dimensions

VSB 25 Blower

VSB 25




Dimensional Data


A	9.92	F	3.68	P	16.25
B	14.6	G	6.60	Y	6.00
C	12.40	H	1.40	Z	11.50
D	7.87	D1	8.62	L	3.80
E	4.12	M	4.20	X ₁	3.25
All dimensions in inches				X ₂	20.75

Wheel Diameter 9.84"



Dimensions for P, X₁, and X₂ will vary depending upon motor selection.



Dimensional Data Roof Unit



A	20.00
B	25.60
C	23.60
D	8.00
E	21.90 x 21.90 (inside)
All dimensions in inches	



R.P.M.	LWA dB (A)
1150	48
1725	58
3450	75